




Figure 1




Digest DNA with  
*Csp6 I* and *Xma I*.



Digest overhangs with  
mung bean exonuclease.




Ligate *Sca I* linkers  
in the presence of *Sca I*.




**Hybridize fragments  
to a biotin labeled probe**


**Figure 2**



**Add blocked streptavidin  
coated magnetic beads.**




**Magnetic capture the magnetic beads,  
biotin labeled probe and associated fragments.**



**Figure 3**

**7. Wash beads and elute DNA.**



**Clone and sequence to characterize captured fragments.**